

new!

# cavity margin insert tips

Preparation of  
crown margins –  
controlled and  
without bleeding



cavity margin insert tips allow a defined preparation border on the gingival margin

The new cavity margin insert tips offer the solution to the difficult field of crown preparation. Their special torpedo shaped contour allows optimal positioning and a precisely defined crown margin. Their oscillating movement is gentle on the gums – no more bleeding – and makes it easier to get the impression.



mectron Srl

Via Loreto, 15/A  
16042 Carasco (Ge)  
ITALIA  
tel +39 0185 35361  
fax +39 0185 351374

www.mectron.com  
mectron@mectron.com

smile. we like it!

# cavity margin insert tips

Order by fax at +39 0185 351347

	<b>CM1</b>	Reference number 02190003	Price (€)* 106,00	Quantity
	Shape Angled shaft with torpedo shaped working tip, Ø 1.8 mm, medium diamond coating (D91)	Treatment Preparation of the crown core and the crown margin		
	<b>CM2</b>	Reference number 02190004	Price (€)* 106,00	Quantity
	Shape Angled shaft with torpedo shaped working tip, Ø 1.8 mm, fine diamond coating (D30)	Treatment Fine finishing of the crown core and the crown margin		
	<b>CM3</b>	Reference number 02190001	Price (€)* 106,00	Quantity
	Shape Angled shaft with torpedo shaped working tip, Ø 1.5 mm, medium diamond coating (D91)	Treatment Preparation of the crown core and the crown margin		
	<b>CM4</b>	Reference number 02190002	Price (€)* 106,00	Quantity
	Shape Angled shaft with torpedo shaped working tip, Ø 1.5 mm, fine diamond coating (D30)	Treatment Fine finishing of the crown core and the crown margin		

\*recommended sale prices plus local taxes

Surgery stamp

Date, signature